




ROMIDAN
DENTAL QUALITY PRODUCTS



- Cable Connection Test →←
- Demo Mode | **Demo**
- Measurement Cable Icon 
- Display of Selected Numerical Values: 2.0, 1.75, 1.5, 1.25, 1.0, 0.75, 0.5, 0.25, 0.0, FO (Foramen Over)
- High precision of apical foramen localization in wet and dry canals
- Audio feedback with volume control
- Fully automatic measurements
- Compact user friendly design
- Automatic turn-off



RomiApex™ A-15


Romidan's Innovative Apex Locators

RomiApex™ A-15

Romidan's Innovative Apex Locators



Front Panel Available in 3 colors



Accessories

cat. RA15-003



1pcs

Measuring Cable

cat. RA15-004



2pcs

File Clip

cat. RA15-006



5pcs

Lip Clip


cat. RA15-007



1pcs

Touch Probe

cat. RA15-009



Cradle

Specifications:

Dimensions: W55 x H92 x T16mm

Display Type: Color LCD

Power Supply:
1.5V, AAA, Alkaline Battery x1

Display Size: 51 x 38 mm

Weight: 100gr.



Distributed by:



Romidan Ltd.

5, Simcha Holzberg St.
Kiryat-Ono 55022, Israel
Tel: +972-3-7360099
Fax: +972-3-7360097
export@romidan.com
www.romidan.com